## GO-ON® Product Information

GB

### Description

GO-ON® is a sterile non-pyrogenic solution of Sodium Hyaluronate.

### Composition

A syringe of GO-ON® contains 2.5 ml of a 1% Sodium Hyaluronate solution as well as sodium chloride, sodium monohydrogen phosphate, sodium dihydrogen phosphate and water for injection.

### Characteristics

Sodium Hyaluronate is a natural high molecular weight polysaccharide composed of a linear chain of disaccharide units made by sodium glucuronate and N-acetylglucosamine. Sodium Hyaluronate is widely distributed throughout the tissues of the human body and is present in high concentrations in the synovial fluid, of which it is the major component. Sodium hyaluronate plays a significant role in modulating the physical and mechanical interaction between adjacent tissues functioning as a tissue lubricant; it acts also as a viscoelastic support maintaining separation between tissues. Different Sodium Hyaluronate preparations may have different molecular weights but they have the same chemical structure. GO-ON is a 1% Sodium Hyaluronate solution obtained from *Streptococcus equi* by fermentation processes and subsequent purification.

### Indications

GO-ON is indicated as a viscoelastic supplementation for synovial fluid in the knee. Viscosupplementation is recommended in patients with normal weight or moderately overweight with symptomatic mild to moderate knee OA who are not sufficiently improved by nonpharmacological treatment and analgesics/ NSAIDs or are unable to take analgesics/ NSAIDs. The actions of the product are lubrication and mechanical support and it is suited for treatment of the symptoms of osteoarthritis.

### Contraindications

GO-ON® must not be used in patients with a history of hypersensitivity to any of the ingredients of the product and in case of inflammatory joint disease such as rheumatoid arthritis or Bechterew disease. As the safety of this product has not been established in children, pregnant or lactating women, the product must not be used in these patients. Recognized contraindications to the use of intra-articular hyaluronic acid injection include active joint infections, bacteremia and overlying skin disease.

### Warning and Precautions

**This product should only be used by a physician. In particular, the physician should possess training as an Orthopedist, Rheumatologist or Physiatrist.**

Recommendations regarding injection techniques include:

- Ensuring the joint does not show clinical signs of acute inflammation (heat, redness, swelling) before injecting GO-ON;
- Removing excessive synovial fluid (due to joint effusion), if present, before injecting;
- Use aseptic technique to avoid contamination and infection in the injection site.

Do not use if the blister packaging has been opened and/or is damaged. The product must not be used after the expiry date indicated on the box. The product must be kept out of the reach of children. This product may worsen local inflammation in patients with osteoarthritis who present an inflammation of the joint. Therefore, its is advisable to administer the product after the local symptom has been alleviated.

Local pain may rarely occur after administration of the product. To avoid pain after the injection, an instruction should be given to keep the local site as still as possible for 2-6 hours after the injection has been made.

The product must be injected strictly into the intra-articular space otherwise it may cause pain. Avoid injection into a blood vessel or surrounding tissues! In the case of patients with a history of hypersensitivity to drugs, administer with caution and keep the patient under clinical monitoring for at least 1 hour after injection.

### Side effects

Some patients may rarely develop a skin rash, such as urticaria and pruritus. In such cases, discontinue the administration and give the patient the appropriate treatment. Pain (mainly transient pain after administration), swelling and effusion, as well as infection, at the injection site have rarely been reported. Patients may rarely complain of redness, heat sensation and heaviness at the injection site. Shock has been very rarely reported. Anaphylactic reactions may occur. If this occurs, discontinue administration, closely monitor the patient and provide adequate treatment. Other reported incidents include reactive arthritis, vastus medialis and skin necrosis, pulmonary embolism, lesions, ruptured meninges when injected into facet joints, cutaneous eruption, deep phlebitis, general fatigue, hypertension, high blood pressure, cardiac arrhythmia, knee infection, popliteal cyst, allergic reaction and surgical intervention. The injection of GO-ON® may be followed by local side effects in the treated joint, e.g. pain,

a sensation of heat, reddening and/or swelling and effusion. These adverse events can be reduced by the application of an ice-bag to the treated joint for 5 to 10 minutes.

### Interactions with Other Agents

There is a known incompatibility between sodium hyaluronate and quaternary ammonium salts such as benzalkonium chloride which can be used as skin disinfectants. GO-ON® must therefore never come into contact with such products, nor with medical or surgical equipment treated with these types of products, since hyaluronic acid can precipitate in their presence. To date there are no reports of incompatibilities (interactions) of GO-ON® with other solutions for intra-articular administration. The patient may benefit from the concurrent use of oral analgesic and anti-inflammatory agents during the first days of treatment.

### Dosage and Administration

GO-ON should be injected into the affected joint a total of 5 times, at weekly intervals. Depending on the severity of the joint disease, the effects of a course of five treatments may last for over 6 months. Treatment cycles may be repeated if required. In the event of an effusion into a joint, it is recommended to carry out a puncture of the effusion, to immobilize the joint, apply an ice bag and/or to administer corticosteroids by intra-articular injection. The treatment with GO-ON® may be repeated 2-3 days later. Remove the syringe from the blister, remove the rubber stopper connection, mount a suitable sterilized needle (e.g. 19 or 21 G) and secure it by a slight turn. Do not over tighten or apply excessive leverage when attaching the needle or removing the needle guard. Remove the air bubble, if any, from the syringe before the injection. GO-ON should be administered on its own and not mixed with any other substances prior to injection.

### Properties and Mode of Action

Synovial fluid, which due to its content of Hyaluronic acid is visco-elastic, occurs in all synovial joints but especially in the large weight bearing joints, where its lubricating and shock absorbing characteristics ensure normal, painless movement. It also supplies nutrients to the articular cartilage. Degenerative joint diseases such as osteoarthritis are associated with a substantial loss of viscosity of the synovial fluid, which impairs its lubricating and shock absorbing functions. This increases the mechanical stress on the joint as well as the loss of articular cartilage to such an extent as to cause pain and loss of function in the affected joints. It has been established that an improvement in the quality of the synovial fluid resulting from an intra-articular injection of Sodium Hyaluronate preparations improves the visco-elastic properties of the synovial fluid. This enhances its lubricant as

well as its shock absorbing actions and reduces the mechanical stress on the joint. This results in alleviation of the pain and an improvement in joint mobility, which after a single treatment cycle of 5 intra-articular injections, last for minimum of 6 months.

### Packaging

GO-ON® is available as a disposable syringe containing a 2.5 ml solution of Sodium Hyaluronate (25 mg). The product is sterilized using steam within blister. Packages contain 1 syringe.

### Shelf-life

Three years

### Storage conditions



Store at a temperature between 2°C and 25°C

### Symbols on the package


Batch number	
Use by	
Use only once	
Attention, see Product Information	
Sterilised by moist heat	
Storage temperature 2°C-25°C	
Manufacturer	
Do not open if the package is damaged	
Keep dry	
Do not resterilize	

### HOW TO USE GO-ON®


Take the syringe out of the container. In case the product is stored in the fridge (see storage conditions), take it out twenty or thirty minutes before its use.


Remove the rubber stopper by turning it

Put on a 19 or 21 G injection needle. Do not over tighten or apply excessive leverage when attaching the needle or removing the needle guard, as this may break the syringe tip.

The syringe is ready to use

\*Before the administration, disinfect the injection site thoroughly. Do not concomitantly use disinfectants containing quaternary ammonium salts.

Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.

Please dispose of GO-ON in accordance with local regulations.

Note: The GO-ON implant is considered safe in the Magnetic Resonance (MR) Environment.

### Manufacturer:

Rottapharm Ltd.,  
Damastown Industrial Park,  
Mulhuddart-Dublin 15  
Ireland  
+353 1 8852 700  
0123

### Distributor for (Country):

Please see carton box

This leaflet was last issued 09/2022

M1045N

4042

<b>ROTTAPHARM Ltd.</b>		Part Number: <b>M1045N</b>	Date: <b>13/10/2022</b>
<b>Proof number:</b> 03	<b>Pharmacode:</b> 4042	<b>Inks:</b> ▼ BLACK	
<b>Job description:</b> Leaflet Go-On Eastern Europe	<b>Barcode format &amp; number:</b> N/A		
<b>Keyline reference:</b> 450 X 210 mm	<b>Braille text:</b> N/A	<b>Template:</b> ---KEYLINE	
<b>Font Used:</b> Arial (8pt)	<b>Artworker:</b> D. Varga		
<b>Name:</b> _____	<b>Approved &amp; proceed</b> <input type="checkbox"/>	<b>Signature:</b> _____	
<b>Date:</b> _____	<b>Amend &amp; re-submit</b> <input type="checkbox"/>		
developed by			

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